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Literatura

1. Initial anatomic investigations of the I-gel airway: a novel supraglottic airway without inflatable cuff (R. M. Levitan and W. C. Kinkle). *Anaesthesia*, 2005, 60, pages 1022-1026.
2. The i-gel supraglottic airway: A potential role for resuscitation? (D Gabbott and R Beringer). *Resuscitation* (2007) 73, 161-164.
3. The i-gel airway for ventilation and rescue intubation. (S. Sharma. R. Rogers. M. Popat). *Anaesthesia* 2007, 62. Correspondence. pages 412-423.
4. Evaluation of four airway training manikins as patient simulators for the insertion of eight types of supraglottic airway devices. (K. M. Jackson¹ and T. M. Cook). *Anaesthesia*, 2007, 62, pages 388-393.
5. The i-gel supraglottic airway and resuscitation - some initial thoughts (J Soar) *Resuscitation* (2007), doi: 10.1016/j.resuscitation. 2007 02.012
6. R. M. Levitan and W. C. Kinkle. Initial anatomic investigations of the I-gel airway: a novel supraglottic airway without inflatable cuff. *Anaesthesia* 2005 Oct;60(10):1022-6.
7. Miller, DM. A proposed classification and scoring system for supraglottic sealing airways: A brief review. *Anesthesia and Analgesia* 2004; 99: 1553-1559.
8. Brimacombe J, Berry A. The incidence of aspiration associated with the laryngeal mask airway-a meta-analysis of published literature. *J Clin Anesth* 1995; 7: 297-30.
9. Brimacombe J, Keller C. Aspiration of gastric contents during use of a ProSeal laryngeal mask airway secondary to unidentified foldover malposition. *Anesthesia Analgesia* 2003; 97: 1192-4 *British Journal of Anaesthesia*
10. Cook TM, Nolan JP, Verghese C, Strube PJ, Lees M, Millar JM, Baskett PJ: Randomized crossover comparison of the Proseal with the Classic laryngeal mask airway in unparalysed anaesthetized patients. *British Journal of Anaesthesia* 2002;88:527-33
11. Brimacombe J, Keller C, Fullekrug B, Agro F, Rosenblatt W, Dierdorf SF, Garcia de Lucas E, Capdevila X, Brimacombe N: A multicenter study comparing the ProSeal with the Classic laryngeal mask airway in anesthetized, nonparalyzed patients. *Anesthesiology* 2002; 96:289-95
12. Brimacombe J, Keller C: The ProSeal laryngeal mask airway: A randomized, crossover study with the standard laryngeal mask airway in paralyzed, anesthetized patients. *Anesthesiology* 2000; 93:104-9
13. Brain AJ, Verghese C, Strube PJ: The LMA 'ProSeal': A laryngeal mask with an oesophageal vent. *British Journal of Anaesthesia* 2000; 84:650-4
14. Keller C, Brimacombe J, Keller K, Morris R. A comparison of four methods for assessing airway sealing pressure with the laryngeal mask airway in adult patients. *British Journal of Anaesthesia* 1999; 82: 286-7
15. Brimacombe J, Keller C, Giampalmo M, Sparr HJ, Berry A. Direct measurement of mucosal pressures exerted by cuff and non-cuff portions of tracheal tubes with different cuff volumes and head and neck positions. *British Journal of Anaesthesia* 1999; 82: 708-11

16. Brimacombe J, Keller C. A comparison of pharyngeal mucosal pressure and airway sealing pressure with the laryngeal mask airway in anesthetized adult patients. *Anesthesia & Analgesia* 1998; 87:1379-82.
17. Keller C, Puehringer F, Brimacombe J. The influence of cuff volume on oropharyngeal leak pressure and fiberoptic position with the laryngeal mask airway. *British Journal of Anaesthesia* 1998;81:186-7.
18. Brimacombe J, Keller C, Morris R, Mecklem D. A comparison of the disposable versus the reusable laryngeal mask airway in paralyzed adult patients. *Anesthesia & Analgesia* 1998;87:921-24.
19. Twigg S, Brown JM, Williams R. Swelling and cyanosis of the tongue associated with use of a laryngeal mask airway. *Anaesthesia & Intensive Care*. 2000; 28:449-50.
20. Stewart A, Lindsay WA. Bilateral hypoglossal nerve injury following the use of the laryngeal mask airway. *Anaesthesia*. 2002; 57:264-5.
21. Drummond GB. Influence of thiopentone on upper airway muscles. *British Journal of Anaesthesia* 1989; 63: 12-21
22. Lowinger D, Benjamin B, Gadd L. Recurrent laryngeal nerve injury caused by a laryngeal mask airway. *Anaesthesia & Intensive Care*. 1999; 27:202-5.
23. Ouellette RG. The effect of nitrous oxide on laryngeal mask cuff pressure. *American Association of Nurse Anesthetists Journal*. 2000; 68:411-4.
24. Agro F, Cataldo R, Carassiti M, Costa F. The seeing stylet: A new device for tracheal intubation. *Resuscitation*. 2000; 44:177-80.
25. Levitan RM, Ochroch AE, Hollander J, et al. Assessment of Airway Visualization: Validation of the Percent of Glottic Opening (POGO) Scale. *Academic Emergency Medicine*. 1998; 5: 919-23.
26. Ochroch AE, Kush S, Stuart S, Hollander JE, Levitan RM. Assessment of laryngeal view in direct laryngoscopy: The percentage of glottic opening (POGO) score compared to Cormack and Lehane grading. *Canadian Journal of Anesthesia*. 1999; 46: 987-90.
27. Cooper RM. Use of a new videolaryngoscope (GlideScope) in the management of a difficult airway. *Canadian Journal of Anesthesia*. 2003; 50:611-3.
28. Brimacombe J, Berry A. A proposed fiber-optic scoring system to standardize the assessment of laryngeal mask airway position. *Anesthesia Analgesia* 1993; 76: 45
29. Keller C, Brimacombe J: Mucosal pressure and oropharyngeal leak pressure with the ProSeal versus the classic laryngeal mask airway. *British Journal of Anaesthesia* 2000; 85:262-6
30. Brimacombe J, Keller C, Boehler M, Puehringer F: Positive pressure ventilation with the ProSeal versus Classic laryngeal mask airway: A-3 randomized, crossover study of healthy female patients. *Anesthesia Analgesia* 2001;93:1351
31. Preoperative fasting for adults to prevent perioperative complications. Brady M, Kinn S, Stuart P. *Cochrane Database Syst Rev*. 2003;(4):CD004423
32. Brimacombe J, *Laryngeal Mask Anaesthesia, Principles and Practice* 2nd edition 2004

All rights reserved. No part of this publication may be reproduced, stored in a retrieval system or transmitted in any form or by any means electronic, mechanical, photocopying, recording or otherwise, without the prior written permission of the publisher.

i-gel is a registered trademark of Intersurgical. The information provided in this document are correct at the time of publication. Intersurgical reserves the right to improve or modify the product without prior notification

Intersurgical Limited
 The Parkway, Warrington, Cheshire, WA1 2RQ
 UK
 Tel: +44 (0)1925 838300
 Fax: +44 (0)1925 838301
 Email: enquiries@intersurgical.com
www.intersurgical.com

Deutschland
 Intersurgical GmbH
info@intersurgical.de

France
 Intersurgical SAS
info@intersurgical.fr

India
 Intersurgical Pvt. Ltd.
info@intersurgical.in

Italy
 Intersurgical S.p.A.
info@intersurgical.it

Japan
 Intersurgical (Japan) Ltd.
info@intersurgical.jp

Spain
 Intersurgical S.L.
info@intersurgical.es

Sweden
 Intersurgical AB
info@intersurgical.se

Switzerland
 Intersurgical AG
info@intersurgical.ch

USA
 Intersurgical Inc.
info@intersurgical.com

USA
 Intersurgical Inc.
info@intersurgical.com

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K-Y Jelly® is a registered trademark of Johnson and Johnson Inc.